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
DÉPARTEMENT INTERNATIONAL
INTERNATIONAL DEPARTMENT

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

4. Cleaning of medical devices	44
4.1 Why should all medical devices be thoroughly cleaned before processing?	44
4.2 Point-of-use preparation of devices for decontamination	44
4.3 Fundamental role of cleaning	45
What is cleaning?	46
Effects of not dismantling and cleaning devices	46
Factors that affect cleaning	46
4.4 Cleaning products	47
Properties associated with ideal cleaning agents	47
Selection of cleaning agents	47
Enzymatic (proteolytic) cleaners	47
Cleaning chemicals (detergents)	47
Preparation of detergent solutions	48
Lubricants	49

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
Source: Health Building Note 13 (HBN13), Department of Health, United Kingdom, 2004

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Cleaning

High Level
Disinfection

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Post-procedural sorting, separating and accounting for contents of sets of devices is most important for several reasons.

- To ensure that all devices and their parts are present prior to cleaning
- To keep reusable devices in designated areas
- To inspect devices for signs of damage and make a note of any repairs necessary
- To allow for tracking and traceability of single devices or sets of devices/instruments, particularly in the case that the infection prevention and control team or operating department need to know on whom the device or devices have been used, e.g. concerns associated with Creutzfeldt-Jakob disease and other risks
- Transport of contaminated devices post-procedure should be sent to a designated decontaminated area as soon as possible and secured in closed, leak and puncture-resistant containers/trolley with tamper-proof locks and security tags
- Biohazard sign should be clearly visible
 - Internationally-recognized symbol or biohazard signage clearly identified in different colours – orange, yellow, red

There are various types of closed container systems available. These include trays, trolleys, impermeable bags, lidded bins and rigid container systems.

- Transport containers should protect both the equipment and the handler from accidental contact with blood or body fluids during transit. These containers may only be reused following cleaning
- Ideally, devices should be transported in moist conditions, without excess liquid

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
ANIOSYME PRIME

Enzymatic complex

Observation of proteasic activity of the cleaning disinfectant ANIOSYME PRIME and put in touch during 15 minutes with the substrate human albumin.


ANIOSYME PRIME

Enzymatic complex

Observation of lipasic activity of the cleaning disinfectant ANIOSYME PRIME put in touch during 15 minutes with the substrate triglycerids


ANIOSYME PRIME

Enzymatic complex

Observation of amylasic activity of the cleaning disinfectant ANIOSYME PRIME put in touch during 15 minutes with the substrate glycogen

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
ANIOSYME PRIME

Microbiological properties

BACTERICIDIE *Bactericidal activity*

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
ANIOSYME PRIME

Disinfectant
detergent
for instrumentation

- Formulation combining 2 patented concepts of efficiency
- Prevents soil from drying
- Antimicrobial efficacy from 15 min
- Dense foam: instruments are covered evenly
- Wide compatibility

INDICATIONS

Pre-treatment of medical and surgical instruments under foam : during transport in non-humid mode or during storage of coiled instruments before manual or machine cleaning. The product maintains an adequate level of humidity to avoid blood and other organic materials drying.

CHARACTERISTICS

- Tri-enzymatic complex: proteases, lipases, amylases, combined with surfactants
- Non-corrosive for materials
- No production of aerosols, avoiding any risk of emanation from the product
- Ready-to-use, ergonomic foam gun shape
- pH of product : 8 +/- 0.5
- Helps the removal of stain when rinsing
- Enzymatic stability proved
- High rinseability

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